

UNIVERSAL SUBMICRON HYBRID COMPOSITE

BRILLIANT EverGlow™

The glow of the art



PERMANENT GLOSS IN NEXT TO NO TIME

BRILLIANT EverGlow represents a new generation of composites from COLTENE. Due to the intelligent composition including special fillers, this Universal Submicron Hybrid Composite offers an exceptionally uniform and smooth surface. In addition, the versatile filling material excels through its extremely durable gloss properties.



AESTHETIC RESTORATION WITH ONLY ONE SHADE

The development of BRILLIANT EverGlow is based on many years of experience in the field of efficient and effective shade management. The adaptive "duo shade" concept is further improved thanks to state of the art filler-technologies. Due to special blend-in properties, the shades integrate extremely harmoniously into the existing dental arch. A natural-looking restoration is readily achieved with the use of a single shade.

« SEVEN UNIVERSAL AND TWO TRANSLUCENT SHADES FOR AN EASY SHADE SELECTION. »

SMOOTH CONSISTENCY FOR GOOD WETTING ABILITY


The smooth BRILLIANT EverGlow composite can be easily applied into cavities of all classes. In addition, the non-slumping material can be modelled without effort until polymerisation. Due to its good wetting properties it offers excellent adhesion to bonded tooth surfaces without sticking to instruments. The abrasion-resistant material is perfectly suited for aesthetic anterior and posterior restorations.

- Quickly and easily polished to high lustre
- Outstanding gloss retention
- Intuitive colour selection
- Smooth consistency



GLOSS RETENTION


After tooth brush abrasion, measured in gloss units



Source: Internal testing method: 6000 cycles, RDA 75
* not a trade mark of COLTENE

COMPRESSIVE STRENGTH

Measured in MPa



Source: Internal test
* not a trade mark of COLTENE

« BRILLIANT PERFORMANCE »



Initial situation of carious 26 molar



Already shiny surface after formation with BRILLIANT EverGlow and subsequent finishing



End result after final polishing

« BETTER THAN MOST COMPOSITES IN THE MARKETPLACE. »

(D.D.S. FRANK J. MILNAR, USA)

Source: Private Communication

ORDER INFORMATION

BRILLIANT EverGlow Trial Kit

Trial Kit Syringes (4 Syringes) **6001 9720**
A2/B2, A3/D3, A3.5/B3, Translucent (Trans) 3 g each

Trial Kit Tips (32 Tips) **6001 9721**
A2/B2, A3/D3, A3.5/B3, Translucent (Trans) 8 x 0.2 g each

BRILLIANT EverGlow System Kit

System Kit Syringes (9 Syringes + Accessory) **6001 9722**
Bleach (BL), A1/B1, A2/B2, A3/D3, A3.5/B3, C2/C3, A4/C4, Translucent (Trans), Bleach Translucent (BL Trans) 3 g each

System Kit Tips (60 Tips + Accessory) **6001 9723**
A1/B1, A2/B2, A3/D3, A3.5/B3, Translucent (Trans) 8 x 0.2 g each
Bleach (BL), C2/C3, A4/C4, Bleach Translucent (BL Trans) 5 x 0.2 g each

BRILLIANT EverGlow Anterior Kit

Anterior Kit Syringes (4 Syringes) **6001 9724**
Bleach (BL), A1/B1, Translucent (Trans), Bleach Translucent (BL Trans) 3 g each

Anterior Kit Tips (32 Tips) **6001 9725**
Bleach (BL), A1/B1, Translucent (Trans), Bleach Translucent (BL Trans) 8 x 0.2 g each



BRILLIANT EverGlow Refills Syringes (3 g each)


| | |
|-------------------------------------|------------------|
| Bleach (BL) | 6001 9697 |
| A1/B1 | 6001 9698 |
| A2/B2 | 6001 9699 |
| A3/D3 | 6001 9700 |
| A3.5/B3 | 6001 9701 |
| C2/C3 | 6001 9702 |
| A4/C4 | 6001 9703 |
| Translucent (Trans) | 6001 9704 |
| Bleach Translucent (BL Trans) | 6001 9705 |

BRILLIANT EverGlow Refills Tips (8 x 0.2 g each)

| | |
|-------------------------------------|------------------|
| Bleach (BL) | 6001 9706 |
| A1/B1 | 6001 9707 |
| A2/B2 | 6001 9708 |
| A3/D3 | 6001 9709 |
| A3.5/B3 | 6001 9710 |
| C2/C3 | 6001 9716 |
| A4/C4 | 6001 9717 |
| Translucent (Trans) | 6001 9718 |
| Bleach Translucent (BL Trans) | 6001 9719 |



© 2015 Coltène/Whaledent AG - www.coltene.com



Coltène/Whaledent AG
Feldwiesenstr. 20
9450 Altstätten / Switzerland
T +41 71 757 53 00
F +41 71 757 53 01
info.ch@coltene.com

